



# DCR Intubation Stent\*



Provided with a molded tip<sup>3</sup> at the insertion end of the stent and a large tapered-shaft punctum plug<sup>2</sup> with a slightly oversized rim molded into the tubing. The metal stylet has a plastic handle<sup>1</sup> for ease of insertion.



- 40.0mm overall length – tubing O.D. 0.79mm tubing I.D. 0.50mm
- 1.0mm plug diameter
- Patency of the stent allows irrigation after insertion

## DESCRIPTION

## REF

DCR Intubation Stent

1087

Supplied: 1 each sterile silicone stent per single pack, preloaded, on a disposable stylet.

\*Developed by Joshua Frankel, MD, Columbia, SC.

