


Graether 2000™

PUPIL EXPANDER SYSTEM

The Graether Pupil Expansion System gently dilates the difficult-to-manage pupil to a natural circular opening of 6.3mm.



Supplied


1 sterile Iris Glide Retractor,
1 sterile silicone ring for
dilation preloaded on
an insertion instrument.

Description

Graether 2000 Pupil Expander

REF

1550



CE Marked

- Gentle Circular Dilation
- 6.3mm Aperture
- Pre-loaded
- Ideal For:
 - Intraoperative Floppy Iris Syndrome
 - Small Pupil
- Instructional video available

